


stroke treatment

corail+

guiding catheter
with balloon



Corail+

Reference	I.D.	O.D.	Long tip length (cm)	Total length (cm)	Ø Inflated balloon 1 mL (mm)	Max balloon volume (mL)	No. of RX markers
CORAIL6F+	1,20 mm (.047")	6F	/	100	11	1,20	1 distal to the balloon
CORAIL8F+	1,89 mm (.074")	8F	/	100	12	1,20	1 distal to the balloon
CORAIL8F+LT30	1,89 mm (.074")	8F	3	110	12	1,20	1 distal to the long tip